



## ***Washing Buffer PBS without detergents***

### **For use in immunohistochemistry**

#### Product description

Composition:	PBS-based buffer with 350 mM NaCl 10x stock solution contains neither Tween nor other detergents	
pH-value:	7.4 ± 0.2	
Preservative:	contains 0.045% ProClin® 300 (v/v)	
Storage:	10x stock solution:	
	shelf life (when stored unopened at -20°C):	1 year - tolerates repeated freezing and thawing cycles
	shelf life (when stored unopened at 2 - 8 °C):	6 months
Use working solution immediately!		

**For research use only, not for diagnostic use**

#### Available

Package sizes: 500 mL      order number 141 500

#### Instructions for use

*Washing Buffer PBS without Tween* is designed especially for immunohistochemistry. It should be used, whenever detergents can make problems in detection.

Crystals of salt can precipitate after storage at 2 - 8 °C or after freezing. Therefore *Washing Buffer* must be warmed up to room temperature and should be mixed thoroughly before preparing the working solution. This leads to dissolving of salt after shaking. Stock solution is diluted 1:10 with salt free water to get the working solution.  
Use working dilution immediately.

For further information please visit [www.candor-bioscience.com](http://www.candor-bioscience.com).

ProClin is a registered trade mark of Rohm and Haas Company.