



## Mouse anti-CD34

Cat. No.: BMS045 (16 ml Ready-to-use)

### Instructions for use

#### Intended use

This antibody is designed for the specific localisation of CD34 in formalin-fixed, paraffin-embedded tissue sections, in frozen sections as well as in cytological specimens.

Anti-CD34 antibody is intended for in vitro diagnostic use.

#### Specifications

<b>Specificity:</b>	CD34
<b>Immunogen:</b>	Endothelial cell membranes obtained as vesicles from human placenta
<b>Clone:</b>	QBEnd/10
<b>Isotype:</b>	Mouse IgG1 kappa
<b>Species reactivity:</b>	Human +, others not tested

#### Summary and Description

CD34 is useful for the immunohistochemical identification of vascular and lymphatic tumours and for the sub-classification of leukaemia.

CD34 is a cell surface protein with a molecular weight of approximately 116 kDa, expressed on hematopoietic stem/progenitor cells and the majority of leukaemia. CD34 is also detected in capillary endothelial cells, where proliferating endothelial cells show higher expression levels than non-proliferating cells.

#### Reagent provided

Mouse monoclonal antibody in buffer with carrier protein and preservative for stabilisation in the following formats:

**Ready-to-use** 16 ml (Cat. No. BMS045)

#### Dilution of primary antibody

None

#### Storage and handling

The antibody should be stored at 2-8°C without further dilution.

If necessary, dilutions of the antibody should be done with a suitable antibody dilution buffer (e.g. ZUC025 from Zytomed Systems). The diluted antibody should be stored at 2-8°C after use. Stability of this working solution depends on various parameters and has to be confirmed by appropriate controls.

The antibody provided is suitable for use until the expiry date indicated on the label, if stored at 2-8°C. Do not use product after the expiry date. Positive and negative controls should be run simultaneously with all specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Zytomed Systems' technical support or your local distributor.

#### Precautions

Use through qualified personnel only.

Wear protective clothing to avoid contact of reagents and specimens with eye, skin and mucous membranes. If reagents or specimens come in contact with sensitive area, wash with large amounts of water.

Microbial contamination of the reagent must be avoided, since otherwise non-specific staining may occur.

ProClin300 is used for stabilisation. A Material safety data sheets (MSDS) is available upon request.

#### Staining procedure

Refer to the following table for conditions specifically recommended for this antibody. Also refer to detection system data sheets for guidance on specific staining protocols or other requirements.

<u>Parameters</u>	<u>Zytomed Systems recommendations</u>
*Pre-treatment	optional: Heat Induced Epitope Retrieval (for example in Citrate buffer pH 6.0 ZUC028)
*Control tissue	Tonsil, placenta
*Working dilution	None
*Incubation time	30 - 60 minutes

## Quality control

The recommended positive control tissues for this antibody are tonsil or a placenta. We recommend carrying out a positive and a negative control with every staining run. Please refer to the instructions of the detection system for guidance on general quality control procedures.

## Troubleshooting

If you observe unusual staining or other deviations from the expected results please read these instructions carefully, refer to the instructions of the detection system for relevant information or contact your local distributor.

## Expected results

This antibody stains positive in the cytoplasmic membrane of haematopoietic stem cells and vascular endothelial cells in formalin-fixed, paraffin-embedded tissue sections. Further details about the expression pattern of p40 can be found in the chapter 'Summary and Description'. Interpretation of the staining results is solely the responsibility of the user. Any experimental result should be confirmed by a medically established diagnostic procedure.

## Limitations of the Procedure

Immunohistochemistry is a complex technique involving both histological and immunological detection methods. Tissue processing and handling prior to immunostaining, for example variations in fixation and embedding or the inherent nature of the tissue can cause inconsistent results (Nadji and Morales, 1983). Endogenous peroxidase, alkaline phosphatase or biotin may cause non-specific staining depending on the detection system used. Tissues containing Hepatitis B Surface Antigen (HBsAg) may give false positive results with HRP (horse radish peroxidase) detection systems (Omata *et al*, 1980). Inadequate counterstaining and mounting can influence the interpretation of the results.

Zytomed Systems warrants that the product will meet all requirements described from its shipping date until the expiry date is reached, if the product is stored and utilised as recommended. No additional guarantees can be given. Under no circumstances shall Zytomed System be liable for any damages arising out of the use of the reagent provided.

## Performance characteristics

Zytomed Systems has conducted studies to evaluate the performance of the antibody for use with a standard detection system. The product has been found to be sensitive and specific to the antigen of interest with minimal or no cross-reactivity.

## Bibliography

Salizzoni M *et al*. Transplantation 76:844-848, 2003  
Fanburg-Smith JC *et al*. Mod Pathol 16:263-271, 2003  
Fina L *et al*. Blood 75:2417-2426, 1990  
Sankey EA *et al*. J Pathol 43:752-757, 1990






Ramani P *et al*. Histopathol 17:237-242, 1990  
Nadji M and Morales AR Ann N.Y. Acad Sci 420:134-9, 1983  
Omata M *et al*. Am J Clin Pathol 73(5): 626-32, 1980

March 03, 2014

Rev: A0314

Doc: DBE\_BMS045

Explanations of the symbols on the product label:

<b>REF</b>	Bestellnummer Catalog Number Reference du catalogue		Verwendbar bis Use By Utiliser jusqu'e		Gebrauchsanweisung beachten Consult Instructions for use Consulter les instructions d'utilisation
<b>LOT</b>	Chargenbezeichnung Batch Code Code du lot		Lagerungstemperatur Temperature Limitation Limites de temp'ature	<b>RUO</b>	Nur f'ur Forschungszwecke For Research Use Only Pour la recherche uniquement
<b>IVD</b>	In vitro Diagnostikum In Vitro Diagnostic Medical Device Dispositif m'edical de diagnostic in vitro		Achtung Warning Attention		Hersteller / Manufacturer / Fabricant Zytomed Systems GmbH • Anhaltinerstra'e 16 14163 Berlin, Germany • Tel: (+49) 30-804 984 990 www.zytomed-systems.com