

## CD11c (Integrin Alpha-X) (EP157)


### Rabbit anti-human CD11c (Integrin Alpha-X) Monoclonal Antibody (Clone EP157)

#### REFERENCES AND PRESENTATIONS<sup>1</sup>

- **ready-to-use (manual or LabVision AutoStainer)**  
MAD-000674QD-3  
MAD-000674QD-7  
MAD-000674QD-12
- **Ready-to-use (MD-Stainer)<sup>2</sup>**  
MAD-000674QD-3/V  
MAD-000674QD/V
- **concentrated**  
MAD-000674Q - 1:50 recommended dilution

#### COMPOSITION

Anti-human CD11c rabbit monoclonal antibody purified from serum and prepared in 10mM PBS, pH 7.4, with 0.2% BSA and 0.09% sodium azide

**INTENDED USE** : Immunohistochemistry (IHC) on paraffin embedded tissues. Not tested on frozen tissues or Western-Blotting.

**CLONE:** EP157<sup>3</sup>

**Ig ISOTYPE:** Rabbit IgG

**IMMUNOGEN:** A synthetic peptide corresponding to residues of human CD11c/ITGAX protein.

**SPECIES REACTIVITY:** In vitro diagnostics in humans. Not tested in other species



**DESCRIPTION AND APPLICATIONS:** CD11c (ITGAX), a member of the leukointegrin family, shares the same beta subunit with other members of the leukocyte adhesion molecule family, which includes CD11a (LFA-1), CD11b (MAC-1) and CD11d (ITGAD), but has a unique alpha chain. CD11c has been shown to play a role in phagocytosis, cell migration, and cytokine production by monocytes/macrophages as well as induction of T cell proliferation by Langerhans cells. CD11c is expressed prominently on the plasma membranes of monocytes, tissue macrophages, NK cells, and most dendritic cells (DCs). A lower level of expression is also observed on neutrophils as a result

of its high level of expression on most DCs. An antibody to CD11c may aid in identification of lesions with histiocytic origin. It may also been used as a marker for hairy cell leukaemia in paraffin embedded tissues.

**IHC POSITIVE CONTROL:** Tonsil or hairy cell leukaemia  
**VISUALIZATION:** Cell membrane and cytoplasm

#### IHC RECOMMENDED PROCEDURE:

- 4µm thick section should be taken on charged slides; dry overnight at 60°C
- Deparaffinise, rehydrate and HIER (heat induced epitope retrieval) – boil tissue in the Pt Module using Vitro S.A EDTA buffer pH8<sup>4</sup> for 20 min at 95°C. Upon completion rinse with 3-5 changes of distilled or deionised water followed by cooling at RT for 20 min
- Endogenous peroxidase block - Blocking for 10 minutes at room temperature using peroxidase solution (ref. MAD-021540Q-125)
- Primary antibody: incubate for 10 minutes [The antibody dilution (when concentrated) and protocol may vary depending on the specimen preparation and specific application. Optimal conditions should be determined by the individual laboratory]
- For detection use Master Polymer Plus Detection System (HRP) (DAB included; ref. MAD-000237QK)
- Counterstaining with haematoxylin and final mounting of the slide

**STORAGE AND STABILITY:**  Stored at 2-8°C. Do not freeze.  Once the packaging has been opened it can be stored until the expiration date of the reagent indicated on the label. If the reagent has been stored under other conditions to those indicated in this document, the user must first check its correct performance taking into account the product warranty is no longer valid.

#### WARNINGS AND PRECAUTIONS:

1. Avoid contact of reagents with eyes and mucous membranes. If reagents come into contact with sensitive areas, wash with copious amounts of water.
2. This product is harmful if swallowed.
3. Consult local or state authorities with regard to recommended method of disposal.
4. Avoid microbial contamination of reagents.

<sup>1</sup> These references are for presentation in vials of Low Density Polyethylene (LDPE) dropper. In case the products are used in automated stainers, a special reference is assigned as follows:

- / L: Cylindrical screw-cap vials (QD-3 / L, QD-7 / L, QD-12 / L).  
- / N: Polygonal screw-cap vials (QD-3 / N, QD-7 / N, QD-12 / N).

For different presentations (references / volumes) please contact the supplier.

<sup>2</sup> For Technical specifications for MD-Stainer, please contact your distributor.

<sup>3</sup> CD11c clone EP157 is manufactured using Epitomics's RabMab® technology under U.S. Patent Nos. 5,675,063 and 7,402,409

<sup>4</sup> Ref: MAD-004072R/D



**Vitro S.A.**

Calle Luís Fuentes Bejarano 60 Ed. Nudo Norte Local 3 41020 Sevilla (Spain)  
Tel: +34 954 933 200. [vitr@vitro.bio](mailto:vitr@vitro.bio) ; [www.vitro.bio](http://www.vitro.bio)



Rev.: 2020-09-21

## SAFETY RECOMMENDATIONS










This product is intended for laboratory professional use only. The product is NOT intended to be used as a drug or for domestic purposes. The current version of the Safety Data Sheet for this product can be downloaded by searching the reference number at [www.vitro.bio](http://www.vitro.bio) or can be requested at [regulatory@vitro.bio](mailto:regulatory@vitro.bio).

## BIBLIOGRAPHY

1. Corbi AL, Larson RS, Kishimoto TK, Springer TA, Morton CC. Chromosomal location of the genes encoding the leukocyte adhesion receptors LFA-1, Mac-1 and p150,95. Identification of a gene cluster involved in cell adhesion. J Exp Med. 1988 May 1;167(5):1597-607
2. Miller ML, Fishleder AJ, Tubbs RR. The expression of CD22 (Leu 14) and CD11c (LeuM5) in chronic lymphoproliferative disorders using two-color flow cytometric analysis. Am J Clin Pathol. 1991 Jul;96(1):100-8
3. Paulli M, Rosso R, Kindl S, Boveri E, Marocolo D, Chioda C, Agostini C, Magrini U, Facchetti F. Immunophenotypic characterization of the cell infiltrate in five cases of sinus histiocytosis with massive lymphadenopathy (Rosai-Dorfman disease). Hum Pathol. 1992 Jun;23(6):647-54
4. Tworek JA, Singleton TP, Schnitzer B, Hsi ED, Ross CW. Flow cytometric and immunohistochemical analysis of small lymphocytic lymphoma, mantle cell lymphoma, and plasmacytoid small lymphocytic lymphoma. Am J Clin Pathol. 1998 Nov;110(5):582-9
5. Vardiman JW, Gilewski TA, Ratain MJ, Bitter MA, Bradlow BA, Golomb HM. Evaluation of Leu-M5 (CD11c) in hairy cell leukemia by the alkaline phosphatase anti-alkaline phosphatase technique. Am J Clin Pathol. 1988 Sep;90(3):250-6
6. Beitnes AC, Ráki M, Brottveit M, Lundin KE, Jahnsen FL, Sollid LM. Rapid accumulation of CD14+CD11c+ dendritic cells in gut mucosa of celiac disease after in vivo gluten challenge. PLoS One. 2012;7(3):e33556
7. Beitnes AC, Ráki M, Lundin KE, Jahnsen J, Sollid LM, Jahnsen FL. Density of CD163+ CD11c+ dendritic cells increases and CD103+ dendritic cells decreases in the coeliac lesion. Scand J Immunol. 2011 Aug;74(2):186-94
8. Ráki M, Tollefsen S, Molberg Ø, Lundin KE, Sollid LM, Jahnsen FL. A unique dendritic cell subset accumulates in the celiac lesion and efficiently activates gluten-reactive T cells. Gastroenterology. 2006 Aug;131(2):428-38.

## LABEL AND BOX SYMBOLS

Explanation of the symbols of the product label and box:

	Expiration date
	Temperature limit
	Manufacturer
	Sufficient content for <n> assays
	Catalog number
	Lot code
	Refer to the instructions of use
	Medical product for <i>in vitro</i> diagnosis.
	Material safety data sheet

