

c-erbB-2/HER2/neu

Catalog # 10047.05 or .10 (0,5 ml or 1,0 ml concentrated) 10047.50 (5 ml ready to use format)

INTENDED USE

- *In Vitro:* This product is intended for qualitative immunohistochemistry on normal and neoplastic, formalin-fixed, paraffin-embedded tissue sections, to be analyzed by light microscopy.
- **Description:** This antibody recognizes a c-erbB-2 protein, which is a receptor tyrosine kinase of the c-erbB family. It is closely related in structure to the epidermal growth factor receptor. C-erbB-2 oncoprotein is detectable in a proportion of breast and other adenocarcinomas, as well as transitional cell carcinomas.
- Expected Staining Pattern: Cell membrane
- **Positive Control:** Breast carcinoma

MATERIALS PROVIDED

- c-erbB-2 specific rabbit monoclonal antibody #10047: purified tissue culture supernatant, with 0.09% sodium azide.
- Antibody Concentration: Not known
- Host: Rabbit
- *Epitope:* Not determined
- Species Reactivity: Human
- Clone: SP3
- Ig Isotype: IgG
- *Immunogen:* Recombinant protein encoding extracellular domain of human c-erbB-2.
- Sterility: This product is not sterile.
- References:
 - 1) Ricardo SA et al. HER2 evaluation using the novel rabbit monoclonal antibody SP3 and CISH in tissue microarrays of invasive breast carcinomas. J Clin Pathol 2007 Sep; 60 (9): 1001-5.

MATERIALS REQUIRED, BUT NOT PROVIDED

- Detection system
- Chromogen/substrate system

- Tissue Section Pretreatment: Staining of formalin-fixed tissue sections requires treating in boiling 10mM citrate buffer, pH 6.0 for 10-20 minutes.
- Recommended Pretreatment Kit: HISTOLS®-Pre-Treatment IHC Kit (cat# 30016), containing the reagents below: HISTOLS®-Citrate Buffer (cat#30010) HISTOLS®-Peroxidase Blocking Solution (cat# 30013) HISTOLS®-BBPS (Background Blocking Protein Solution) (cat# 30012)
- **Dilution of Concentrated Antibody:** 1:100 in antibody diluent.
- **Primary Antibody Incubation Time:** 60 minutes at RT.
- Recommended Detection System: HISTOLS[®]-MR (anti-mouse & rabbit Ig) (cat# 30011)
- *Visualization:* To detect antibody binding sites, follow the instructions provided with the chromogen/substrate system.

STORAGE AND STABILITY

This product contains sodium azide and is stable when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.



c-erbB-2, breast carcinoma

METHODS AND PROCEDURES

