


PolyDetector ALK Scarlet Substrate-Chromogen



Inset: IHC of CD20 with ALK Scarlet on tonsil tissue

Intended Use

For In Vitro Diagnostic Use.

Presentation

PolyDetector ALK Scarlet Substrate-Chromogen kit contains an ALK Scarlet chromogen solution, and an ALK Scarlet Buffer-Substrate solution. All the components are buffered, contain stabilizers and a non-azide antimicrobial.

Summary and Explanation

PolyDetector ALK Scarlet Substrate-Chromogen is suitable for use in Immunohistochemical, Immunocytochemical and in situ hybridization procedures using Alkaline Phosphatase Detection Systems. PolyDetector ALK Scarlet Substrate-Chromogen allows for the demonstration of cell antigens or nucleic acids in paraffin-embedded tissues, cryostat sections, cytospins, and cell preparations.

The substrate-chromogen is the final step in the detection portion; it enables the antibody-antigen-enzyme (or nucleic acid-probe-enzyme) complex to be viewed under the light microscope. This occurs when ALK Scarlet, in the presence of Alkaline Phosphatase, gets deposited at the site of the target antigen or nucleic acid, producing a red color that is partially soluble in organic solvents and therefore can be mounted with permanent medium when using the Bio SB ChromoProtector (BSB 0151 – BSB 0156), Bio SB Fast Chromoprotector (BSB 0327-15-BSB 0327-1000), or using fast xylene/alcohol dehydration protocol.

Catalog No.	Buffer-Substrate	Chromogen
BSB 0138	15 ml	1 mL
BSB 0139	50 ml	3 mL
BSB 0140	100 ml	6 mL
BSB 0141	200 ml	12 mL
BSB 0142	1000 ml	50 ml

Storage Store at 2-8° C

Stability

This product is stable up to the expiration date on the product label.

Do not use after expiration date listed on the package label. Temperature fluctuations should be avoided. Store appropriately when not in use.

Adhere to all local laws when disposing of this product.

Preparation of Working Solution:

To prepare a working PolyDetector ALK Scarlet Substrate-Chromogen solution, add 1 drop of the ALK Scarlet Chromogen to 1 mL of the ALK Scarlet Buffer-Substrate and mix. Use this working solution within 5-10 minutes of preparation.

Alk Scarlet Buffer-Substrate	1 ml	2 ml	3 ml
ALK Scarlet Chromogen	1 drop	2 drops	3 drops

Recommended Immunohistochemical Protocol

Step	ImmunoDetector AP
Peroxidase/AP Blocker	5 min.
Primary Antibody	30-60 min.
1st Detection Step	10 min.
2nd Detection Step	10 min.
Substrate-Chromogen	5-10 min.
Counterstain / Coverslip	Varies

Step	PolyDetector AP
Peroxidase/AP Blocker	5 min.
Primary Antibody	30-60 min.
1st Detection Step	30-60 min.
2nd Detection Step	N/A
Substrate-Chromogen	5-10 min.
Counterstain / Coverslip	Varies

Precautions

1. For professional users only. Results should be interpreted by a qualified medical professional.
3. Ensure proper handling procedures are used with reagent. Minimize microbial contamination of reagents.
2. Always wear personal protective equipment such as laboratory coats, goggles and gloves when handling reagents.
3. Dispose of unused solution with copious amounts of water.
4. Do not ingest reagent. If reagent is ingested, seek medical advice immediately.
5. Avoid contact with eyes. If contact occurs, flush with large quantities of water.
6. Follow safety precautions of the heating device used for epitope retrieval (TintoRetriever Pressure Cooker or similar).
7. For additional safety information refer to Safety Data Sheet for this product.
8. For complete recommendations for handling biological specimens, please refer to the CDC document, "Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories" (see References in this document)

Mounting Protocol

a. Alcohol/Xylene Protocol (Fast)

1. After the histological, immunohistochemical or in situ hybridization staining procedure is completed, rinse slides in deionized water.
2. Dip the slides in alcohol 30%, 70%, and 100% for 3-5 seconds, then dip for 3-5 seconds in 3 xylenes.
3. Add an organic Permanent Mounting medium such as XyGreen PermaMounter (BSB 0169-0174), PermaMounter (BSB 0094-0097) or similar permanent mounting media.
4. Apply cover slip and air dry before microscopic observation.

b. ChromoProtector Protocol

1. After the histological, immunohistochemical or in situ hybridization staining procedure is completed, rinse slides in deionized water. Do not incubate tissue or cell specimens in solvents such as alcohol, toluene, or xylene.
2. Lay wet slides horizontally and apply sufficient drops of ChromoProtector (BSB 0151 – BSB 0156) to completely cover the tissue. Carefully spread ChromoProtector if needed, but avoid contacting the tissue.
3. Incubate slides for ten minutes at 60 °C to allow ChromoProtector to penetrate tissues.
4. Remove excess ChromoProtector by placing slides vertically over an absorbent material and let excess drain off into absorbent material. Do not rinse slides.
5. Allow slides to COMPLETELY air dry.

NOTE: The ChromoProtector will protect tissue from drying artifacts during the air-drying process.

6. Add an organic Permanent Mounting medium such as XyGreen PermaMounter (BSB 0169-0174), PermaMounter (BSB 0094-0097) or similar permanent mounting media.
7. Apply cover slip and air dry before microscopic observation.

c. Fast ChromoProtector Protocol

1. Apply a few drops of Fast ChromoProtector onto the tissue. Tilt the slide and let the excess drain off the slide. Let the slide dry completely.
2. Apply an organic permanent mounting medium and cover slip.









Product Limitations

Due to inherent variability present in immunohistochemical procedures (including fixation time of tissues, dilution factor of antibody, retrieval method utilized and incubation time), optimal performance should be established through the use of positive and negative controls. Results should be interpreted by a qualified medical professional.

References

1. U.S. Department of Health and Human Services: Centers for Disease Control and Prevention. Guidelines for Safe Work Practices in Human and Animal Medical Diagnostic Laboratories. Supplement / Vol. 61, January 6, 2012.

Symbol Key / Légende des symboles/Erläuterung der Symbole

 QAdvis EAR AB Ideon Science Park Scheelevägen 17 SE-223 70 Lund, Sweden	 Storage Temperature Limites de température Zulässiger Temperaturbereich	 Manufacturer Fabricant Hersteller	 Catalog Number Référence du catalogue Bestellnummer
 In Vitro Diagnostic Medical Device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	 Read Instructions for Use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	 Expiration Date Utiliser jusque Verwendbar bis	 Lot Number Code du lot Chargenbezeichnung