

Cat No	Item Name & Discription	Package	Stock temperature	Tranport temperature to distributors	Transport temperature to end-customers	Life time
<b>K-12031-C20</b>	WesternBright ECL-HRP Substrate for X-ray Film-sample	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12045-C20</b>	WesternBright ECL- HRP Substrate for X-ray Film-20ml trial kit- contains substrate for 200cm2 membrane	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12045-D20</b>	WesternBright ECL- HRP Substrate for X-ray Film- 200ml kit-contains substrate for 2000cm2	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12045-D50</b>	WesternBright ECL- HRP Substrate for X-ray Flm- 500ml kit-contains substrate for 5000cm2	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12032-C20</b>	WesternBright Quantum- HRP Substrate for CCD- Sample	1kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12042-C20</b>	WesternBright Quantum- HRP Substrate for CCD- 20ml trial kit- contains substrate for 200cm2 membrane	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12042-D10</b>	WesternBright Quantum- HRP Substrate for CCD , 100ml kit-contains substrate for 1000cm2 membrane	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12042-D20</b>	WesternBright Quantum- HRP Substrate for CCD , 200ml kit-contains substrate for 2000cm2 membrane	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12033-C20</b>	WesternBright Sirius- femotogram HRP Substrate , Sample	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12043-C20</b>	WesternBright Sirius- femotogram HRP Substrate , 20ml kit-contains substrate for 200cm2 membrane	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12043-D10</b>	WesternBright Sirius- femtogram HRP Substrate , 100ml kit-contains substrate for 1000cm2 membrane	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>K-12043-D20</b>	WesternBright Sirius- femtogram HRP Substrate , 200ml kit-contains substrate for 2000cm2 membrane	1 kit	RT	RT	RT	RT : 18 MONTHS, 4°C : 2 YEARS
<b>R-03024-D50</b>	AdvanWash- 500ml, Immunoblot Washing Solution 10x concentrate	1 Each	4C	RT	RT	4C/18MONTHS
<b>K-12021-010</b>	WesternBright MFC- Fluorescent Western Blot Kit, kit contains material for 10 Western Blots	1 Kit	4C	RT and 4C	4C	4C/18MONTHS
<b>K-08001-001</b>	Fluorescent Western Calibration Blot, 3 protein fluorescent (Green, Red, Blue channel) WB for imaging instrument calberation	1 Each	RT	RT	RT	3 MONTHS
<b>L-07001-010</b>	Background Quenching Sheets, Proprietary material for Fluorescent gels and WB to reduce auto fluorescence	10 Each	RT	RT	RT	24 MONTHS
<b>R-03023-D20</b>	AdvanBlock-PF- 200ml, Protein-free Blocking Solution (5x concentrate) for Fluorescent Immunoblots	1 Each	4C	RT	RT	18 MONTHS

<b>L-08001-010</b>	Low Fluorescence Western Membrane (PVDF), Pre-cut PVDF membrane (9x7 cm- mini vertical gels)	10 Each	RT	RT	RT	24 MONTHS
<b>L-08002-010</b>	Nitrocellulose Transfer Membrae- 0.45 um, Pre-cut membrane (8x10 cm- mini vertical gels- ideal for larger proteins)	10 Each	RT	RT	RT	18 MONTHS
<b>L-08003-010</b>	Nitrocellulose Transfer Membrae- 0.22 um, Pre-cut membrane (8x10 cm- mini vertical gels- ideal for smaller proteins)	10 Each	RT	RT	RT	18 MONTHS
<b>R-01038-020</b>	Avant Buffer Pouches- PBS, each packet makes 500ml of 1x solution	20/PACK	RT	RT	RT	24 MONTHS
<b>R-01039-020</b>	Avant Buffer Pouches-TBS, each packet makes 500ml of 1x solution	20/PACK	RT	RT	RT	24 MONTHS
<b>R-01037-020</b>	Avant Buffer Pouches-Tris-Glycine PAGE running Buffer, each packet makes 500ml of 1x solution	20/PACK	RT	RT	RT	24 MONTHS
<b>R-01036-020</b>	Avant Buffer Pouches-SDS Tris-Glycine PAGE running Buffer, each packet makes 500ml of 1x solution	20/PACK	RT	RT	RT	24 MONTHS
<b>K-01001-010</b>	RapidClean Protein Removal kit-10 rxns, slurry blue beads and spin filter	1 kit	4C (LIGHT PROTECTED)	RT	RT	24 MONTHS
<b>K-01001-025</b>	RapidClean Protein Removal kit-25 rxns, slurry blue beads and spin filter	1 kit	4C (LIGHT PROTECTED)	RT	RT	24 MONTHS
<b>R-14011-250</b>	RapidClean- 250ul, Slurry blue beads add 1/10 sample volume	1 Each	4C (LIGHT PROTECTED)	RT	RT	24 MONTHS
<b>R-14011-B10</b>	RapidClean- 1ml, Slurry blue beads add 1/10 sample volume	1 Each	4C (LIGHT PROTECTED)	RT	RT	24 MONTHS
<b>K-02101-010</b>	Afyon SDS-PAGE Sample Presperation Kit, Kit contains reagents and columns for 10 protein clean-up	1 kit	4C (LIGHT PROTECTED)	RT	RT	24 MONTHS
<b>K-02101-025</b>	Afyon SDS-PAGE Sample Preparation kit, Kit contains reagents and columns for 25 protein clean-up	1 kit	4C (LIGHT PROTECTED)	RT	RT	24 MONTHS
<b>R-05051-050</b>	APC-Goat-Anti Rrabbit IgG Antibody- 50ul, APC conjugated secodary fluorescent antibody	1 Each	4C	4C	4C	24 MONTHS
<b>R-05051-250</b>	APC-Goat-Anti-Rabbit IgG Antibody- 250ul, APC conjugated secodary fluorescent antibody	1 Each	4C	4C	4C	24 MONTHS
<b>R-05052-050</b>	RPE-Goat-Anti-Mouse IgG Antibody- 50ul, RPE conjugated secodary fluorescent antibody	1 Each	4C	4C	4C	24 MONTHS
<b>R-05052-250</b>	RPE-Goat-Anti-Mouse IgG Antibody- 250ul, RPE conjugated secodary fluorescent antibody	1 Each	4C	4C	4C	24 MONTHS
<b>R-05011-050</b>	AdvanFluor APC-Streptavidin Conjugate- 50ug, APC conjugated streptavidin for protein detection	1 Each	4C	4C	4C	24 MONTHS
<b>R-05011-200</b>	AdvanFluor APC-Streptavidin Conjugate- 200ug, APC conjugated streptavidin for protein detection	1 Each	4C	4C	4C	24 MONTHS
<b>R-05012-050</b>	AdvanFluor BPE-Streptavidin Conjugate- 50ug, BPEconjugated streptavidin for protein detection	1 Each	4C	4C	4C	24 MONTHS

<b>R-05012-200</b>	AdvanFluor BPE-Streptavidin Conjugate- 200ug, BPE conjugated streptavidin for protein detection	1 Each	4C	4C	4C	24 MONTHS
<b>R-05072-500</b>	Goat-anti-rabbit HRP secondary antibody- 500ul, HRP conjugated secondary antibody for chemic Westerns	1 Each	-20C	-20C or 4C	-20C or 4C	24 MONTHS
<b>R-05071-500</b>	Goat-anti-mouse HRP secondary antibody- 500ul, HRP conjugated secondary antibody for chemic Westerns	1 Each	-20C	-20C or 4C	-20C or 4C	24 MONTHS
<b>L-07014-100</b>	LucentBlue X-ray Film- 8x10 cm, optimized blue x-ray film for chemiluminescence	100/PACK	RT	RT	RT	18 MONTHS
<b>L-07013-100</b>	LucentBlue X-ray Film- 5x7 cm, optimized blue x-ray film for chemiluminescence	100/PACK	RT	RT	RT	18 MONTHS
<b>L-07031-005</b>	Incubation Tray- 7x5 cm, for staining and washing gels and membranes	5/PACK	RT	RT	RT	36 MONTHS
<b>L-07032-005</b>	Incubation Tray- 9x6 cm, optimized blue x-ray film for chemiluminescence	5/PACK	RT	RT	RT	36 MONTHS
<b>L-07033-005</b>	Incubation Tray- 11x9 cm- , optimized blue x-ray film for chemiluminescence	5/PACK	RT	RT	RT	36 MONTHS
<b>K-11072-B50</b>	Advanstain Scarlet kit-5ml (20 mini gels)- dye and buffer included in kit for Fluorescent Total Protein Staining	1 Each	-20C	-20C	-20C (can ship 4C only overnight)	18 MONTHS
<b>K-11072-C25</b>	Advanstain Scarlet kit-25ml (100 mini gels)- dye and buffer included in kit for Fluorescent Total Protein Staining	1 Each	-20C	-20C	-20C (can ship 4C only overnight)	18 MONTHS
<b>R-03018-B10</b>	Non-reducing protein sample loading buffer (2X) ready to use laemmli loading buffer	1 Each	RT	RT	RT	18 MONTHS
<b>R-03018-B50</b>	Non-reducing protein sample loading buffer (2X) ready to use laemmli loading buffer	1 Each	RT	RT	RT	18 MONTHS
<b>R-03019-B10</b>	Reducing protein sample loading buffer (2X) ready to use laemmli loading buffer	1 Each	RT	RT	RT	18 MONTHS
<b>R-03019-B50</b>	Reducing protein sample loading buffer (2X) ready to use laemmli loading buffer	1 Each	RT	RT	RT	18 MONTHS
<b>R-03021-D50</b>	Advanstain Ponceau (sufficient for staining 50blots 9x7cm) Visible reversible total protein stain for proteins transferred	1 Each	RT	RT	RT	18 MONTHS
<b>K-16025-C10</b>	ELISABright Substrate- trial sample Chemiluminescent substrate optimized for ELISA( sufficient for 1x96 well plate)	1 Each	RT	RT	RT	18 MONTHS
<b>K-16025-D10</b>	ELISABright Substrate- trial sample Chemiluminescent substrate optimized for ELISA( sufficient for 1x96 well plate)	1 Each	RT	RT	RT	18 MONTHS
<b>K-16025-D25</b>	ELISABright Substrate- trial sample Chemiluminescent substrate optimized for ELISA( sufficient for 1x96 well plate)	1 Each	RT	RT	RT	18 MONTHS
<b>R-07055-001</b>	ChemiPen-Proprietary "ink" for annotating western blots	1 Each	RT	RT	RT	18 MONTHS