


## Product Specifications



Product:	Regeneration Buffer
Description:	Regeneration Buffer
Part Number:	BU-124-FP
Concentration:	n/a
Amount:	50 ml
Form:	Liquid
Purification:	Affinity purified using solid phase biotin
Buffer:	Not Applicable
Storage:	2-8 °C, DO NOT FREEZE
Applications:	For antibody purification
Disclaimer:	For <i>in vitro</i> Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. <b>Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.</b>