

5-methylcytidine/m5C Rabbit pAb

Catalog No.: A18870

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

DB

Cross-Reactivity

Human, Mouse, Rat, Other (Wide Range)

Background

Recommended Dilutions

DB 1:500 - 1:2000

Immunogen Information

Gene ID

Swiss Prot

Immunogen

5mC

Synonyms

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

Validation Data



The m5C rabbit polyclonal antibody (A18870) are tested in Dot Blot against 5-methylcytidine(m5C) and 5-methylcytosine(5mC).