

DLK1 Rabbit pAb

Catalog No.: A6578

Basic Information

Observed MW

45kDa

Calculated MW

33kDa/41kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Background

This gene encodes a transmembrane protein that contains multiple epidermal growth factor repeats that functions as a regulator of cell growth. The encoded protein is involved in the differentiation of several cell types including adipocytes. This gene is located in a region of chromosome 14 frequently showing unparental disomy, and is imprinted and expressed from the paternal allele. A single nucleotide variant in this gene is associated with child and adolescent obesity and shows polar overdominance, where heterozygotes carrying an active paternal allele express the phenotype, while mutant homozygotes are normal.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

8788

Swiss Prot

P80370

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 83-383 of human DLK1 (NP_003827.3).

Synonyms

DLK1;DLK;DLK-1;Delta1;FA1;PREF1;Pref-1;ZOG;pG2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of Mouse heart, using DLK1 Rabbit pAb (A6578) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 180s.