

# ABI3BP Rabbit pAb

Catalog No.: A14581

## Basic Information

### Observed MW

118kDa

### Calculated MW

53kDa/72kDa/86kDa/118kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

25890

### Swiss Prot

Q7Z7G0

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 29-280 of human ABI3BP (NP\_056244.2).

### Synonyms

ABI3BP;NESHBP;TARSH

## Contact

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of extracts of U-87MG cells, using ABI3BP antibody (A14581) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25ug per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 20s.