

TET3 Rabbit pAb

Catalog No.: A13453

Basic Information

Observed MW

Calculated MW

77kDa/166kDa/179kDa

Category

Primary antibody

Applications

IHC, IF/ICC

Cross-Reactivity

Human

Background

Recommended Dilutions

IHC	1:50 - 1:200
IF/ICC	1:50 - 1:200

Immunogen Information

Gene ID	Swiss Prot
200424	O43151

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1556-1795 of human TET3 (NP_001274420.1).

Synonyms

TET3;hCG_40738;TET3

Contact

 | www.abclonal.com

Product Information


Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Validation Data



Immunohistochemistry of paraffin-embedded human colon carcinoma using TET3 antibody (A13453) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human appendix using TET3 antibody (A13453) at dilution of 1:200 (40x lens).



Immunofluorescence analysis of A549 cells using TET3 antibody (A13453). Blue: DAPI for nuclear staining.