

ZZZ3 Rabbit pAb

Catalog No.: A15809

Basic Information

Observed MW

170kDa

Calculated MW

47kDa/76kDa/101kDa/102kDa

Category

Primary antibody

Applications

WB,IF/ICC

Cross-Reactivity

Human, Rat

Background

Component of the ATAC complex, a complex with histone acetyltransferase activity on histones H3 and H4.

Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:100

Immunogen Information

Gene ID	Swiss Prot
26009	Q8IYH5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ZZZ3 (Q8IYH5).

Synonyms

ZZZ3;ATAC1

Contact

 | www.abclonal.com

Product Information


Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using ZZZ3 antibody (A15809) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.