

Phospho-Histone H3-T3 Rabbit pAb

Catalog No.: AP0846

Basic Information

Observed MW

17KDa

Calculated MW

15kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse, Rat, Other (Wide Range)

Background

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6p22-p21.3.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

8350

Swiss Prot

P68431

Immunogen

A synthetic phosphorylated peptide around T3 of human HIST1H3A (NP_003520.1).

Synonyms

H3/A;H3C2;H3C3;H3C4;H3C6;H3C7;H3C8;H3FA;H3C10;H3C11;H3C12;HIST1H3A

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using Phospho-Histone H3-T3 pAb (AP0846) at 1:1000 dilution. C6 and HeLa cells were treated by Nocodazole (50 ng/mL) at 37°C for 20 hours. C6 and NIH/3T3 cells were treated by CIP (20uL/400ul) at 37°C for 1 hour. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 1s.



Western blot analysis of extracts of various cell lines, using Phospho-Histone H3-T3 pAb (AP0846) at 1:1000 dilution or Histone H3 antibody (A2348). Both 293T cells and Raw264.7 cells were treated by Nocodazole (100 ng/mL) at 37°C for 16 hours. C6 cells were treated by Nocodazole (50 ng/mL) at 37°C for 20 hours. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit (RM00020). Exposure time: 1s.