

anti- AARSD1 antibody

Product Information

Catalog No.: FNab00023

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Functions in trans to edit the amino acid moiety from incorrectly charged tRNA(Ala).

Immunogen information

Immunogen: alanyl-tRNAsynthetase domain containing 1

Synonyms: AARSD1
Observed MW: 45 kDa
UniprotID: Q9BTE6

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00023(AARSD1 antibody) at dilution of 1:500

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>