

anti- ABCB7 antibody

Product Information

Catalog No.:	FNab00033
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Could be involved in the transport of heme from the mitochondria to the cytosol. Plays a central role in the maturation of cytosolic iron-sulfur(Fe/S) cluster-containing proteins.

Immunogen information

Immunogen:	ATP-binding cassette, sub-family B(MDR/TAP), member 7
Synonyms:	ABC7, ABCB7, ASAT, Atm1p, EST140535
Observed MW:	83 kDa
Uniprot ID :	O75027

Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IP		
Recommended dilution: WB: 1:500-1:1000; IP: 1:500-1:1000			
Image:			


Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com







mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00033(ABCB7 antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com