

anti- TIMP-1 antibody

Product Information

Catalog No.:	FNab08706
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	TIMP metalloproteinase inhibitor 1
Synonyms:	CLGI, Collagenase inhibitor, EPA, EPO, HCI, Metalloproteinase inhibitor 1, TIMP, TIMP 1, TIMP1
Observed MW:	28 kDa
Uniprot ID :	P01033

Application

Reactivity:	Human, Mouse
Tested Application:	WB, IHC, IF, ELISA
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889


www.fn-test.com



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using FNab08706(TIMP1 Antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells using FNab08706 (TIMP1 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



Various lysates were subjected to SDS PAGE followed by western blot with FNab08706(TIMP1 Antibody) at dilution of 1:800