

## Recombinant Human CD80/B7-1 Protein

Catalog No. PKSH031476

### Description

<b>Synonyms</b>	CD80; Activation B7-1 antigen; B7; BB1; CD28LG1; CD28LGB7-1 antigen; T-lymphocyte activation antigen CD80;B7-1;B7.1;CD28LG;LAB7
<b>Species</b>	Human
<b>Expression_host</b>	HEK293 Cells
<b>Sequence</b>	Met 1-Asn 242
<b>Accession</b>	NP_005182.1
<b>Mol_Mass</b>	25.5 kDa
<b>Tag</b>	None
<b>Bio_activity</b>	Testing in progress

### Properties

<b>Purity</b>	>85 % as determined by reducing SDS-PAGE.
<b>Endotoxin</b>	<1.0 EU per µg of the protein as determined by the LAL method.
<b>Storage</b>	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
<b>Shipping</b>	This product is provided as lyophilized powder which is shipped with ice packs.
<b>Formulation</b>	Lyophilized from sterile PBS; pH 7.4
<b>Reconstitution</b>	Please refer to the printed manual for detailed information.

### Background

The B-lymphocyte activation antigen B7-1 (referred to as B7); also known as CD80; is a member of cell surface immunoglobulin superfamily and is expressed on the surface of antigen-presenting cells including activated B cells; macrophages and dendritic cells. As costimulatory ligands; B7-1 which exists predominantly as dimer and the related protein B7-2; interact with the costimulatory receptors CD28 and cytotoxic T lymphocyte-associated antigen 4 (CTLA-4) expressed on T cells; and thus constitute one of the dominant pathways that regulate T cell activation and tolerance; cytokine production; and the generation of CTL. The B7/CD28/CTLA4 pathway has the ability to both positively and negatively regulate immune responses. CD80 is thus regarded as promising therapeutic targets for autoimmune diseases and various carcinomas.

## SDS-PAGE

