

Recombinant Human Proteoglycan 3/PRG3 Protein (His Tag)

Catalog No. PKSH032976

Description

Synonyms	Proteoglycan 3; Eosinophil Major Basic Protein Homolog; Prepro-Major Basic Protein Homolog; Prepro-MBPH; PRG3; MBPH
Species	Human
Expression_host	Human Cells
Sequence	Leu18-Phe225
Accession	Q9Y2Y8
Mol_Mass	24.6 kDa
AP_Mol_Mass	27 kDa
Tag	C-6His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.2.
Reconstitution	Please refer to the printed manual for detailed information.

Background

Proteoglycan 3, also known as Eosinophil major basic protein homolog, Prepro-major basic protein homolog, PRG3 and MBPH, contains one C-type lectin domain. Proteoglycans are a major component of the animal extracellular matrix. PRG3 localizes to the eosinophil secondary granule and is expressed in bone marrow, not detected in placenta. PRG3 has similar cytotoxic and cyto stimulatory activities to PRG2/MBP. In vitro, PRG3 can stimulate neutrophil superoxide production and IL8 release, histamine and leukotriene C4 release from basophils.

SDS-PAGE

