

Recombinant Rat Tie2/TEK Protein (Fc Tag)(Active)

Catalog No. PKSR030377

Description

Synonyms	Tie-2;Tie2
Species	Rat
Expression_host	HEK293 Cells
Sequence	Met4-Leu743
Accession	NP_001099207.1
Mol_Mass	108 kDa
AP_Mol_Mass	120-130 kDa
Tag	C-Fc
Bio_activity	1. Measured by its binding ability in a functional ELISA.2. Immobilized S1h-3C-mANGPT2 at 10 µg/mL (100 µL/well) can bind ratTEK-Fc.The EC50 of ratTEK-Fc is 0.26-0.62µg/mL.

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4
Reconstitution	Please refer to the printed manual for detailed information.

Background

TEK, or TIE-2, is an endothelial cell-specific receptor tyrosine kinase (RTK) that is known as a functioning molecule of vascular endothelial cells. TEK comprises a subfamily of RTK with TIE, and these two receptors play critical roles in vascular maturation, maintenance of integrity and remodeling. Targeted mutagenesis of both Tek and its agonistic ligand, Angiopoietin-1, result in embryonic lethality, demonstrating that the signal transduction pathways mediated by this receptor are crucial for normal embryonic development. TEK signaling is indispensable for the development of the embryonic vasculature and suggests that TEK signaling may also be required for the development of the tumor vasculature.

SDS-PAGE

