

Cat. No.	1-LI016-07
Lot No.	
Quantity	1 ml

**Antibody to
Lin 28 homolog/LIN28**

Host:	Rabbit
Immunogen:	Synthetic peptide of human Lin28
Presentation:	Purified antibody diluted in TRIS, pH 7.4, containing < 0.1 % sodium azide. Briefly centrifuge small volumes of antibody to remove any liquid from the tube's lid.
Assay system:	IHC(p*,f) *reported as working on paraffin, ELISA, WB
Fixation:	1)NOTOXhisto 2)Formalin
Treatment:	To be determined in user's own system.
Titre:	IHC(p,f) 1:10-1:50, ELISA 1:50-1:100, WB 1:10-1:50; Titre should be considered as a recommendation and depends on incubation time and conditions (heat or RT). Optimum titration to be determined in user's own system. For IHC dilute primary antibody with antibody diluent, Cat# 400100295, and incubate for at least 60 min.
Specificity:	LIN28 is a RNA-binding protein acting as a 'translational enhancer' which drives specific mRNAs to polysomes and thus increases the efficiency of protein synthesis. It controls the timing of events during embryonic development and is readily expressed in embryos, embryonic stem cells and embryonal carcinoma cells. The presence of LIN-28 persists in some adult tissues including cardiac and skeletal muscle.
Reactivity:	Human, Mouse, Rat, others not tested.
Cell. Local.:	Nuclear, Cytoplasmic
Control:	Testis tissue
Storage:	Aliquot and store at - 20 °C. Do not use after expiration date indicated on the vial.

Errors excepted. This data sheet is a general information. The product attribute can diversify with changing Lot No as well as variations in tissue selection, tissue processing, antigen retrieval and detection systems. We do not take responsibility for any possible damage including personal injury, time or effort on economic loss caused by this product. Our warranty is limited to the price paid for the product. The product may only be used by authorized and skilled personnel. Non-application as prescribed in this data sheet leads to loss of all liability.

FOR IN VITRO USE, NOT FOR DIAGNOSTIC USE

Manufacturer:

BIOCYC GmbH & Co. KG, Im Biotechnologiepark TGZ I, 14943 Luckenwalde; certified by

ISO 13485:2012 and AC: 2012, ISO 9001:2008; Registry number: 60018916, 60018917

Distributed by:



Schichauweg 16, 12307 Berlin, Germany
Tel.: ++49 (030) 765 925-0, Fax.: ++49 (030) 765 925-55
E-Mail: info@quartett.com Internet: www.quartett.com