

Recombinant Human CCL5/RANTES Protein

Catalog No.: RP00321 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	6352	P13501

Tags

No tag

Synonyms

D17S136E; eoCP; RANTES; SCYA5; SIS-delta; SISd; TCP228

Product Information

Source	Purification
<i>E. coli</i>	> 95% by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μm filtered solution of 20 mM PB, 150 mM NaCl, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

Background

Chemokines form a superfamily of secreted proteins involved in immunoregulatory and inflammatory processes. The superfamily is divided into four subfamilies based on the arrangement of the N-terminal cysteine residues of the mature peptide. This chemokine, a member of the CC subfamily, functions as a chemoattractant for blood monocytes, memory T helper cells and eosinophils. It causes the release of histamine from basophils and activates eosinophils. This cytokine is one of the major HIV-suppressive factors produced by CD8+ cells. It functions as one of the natural ligands for the chemokine receptor chemokine (C-C motif) receptor 5 (CCR5), and it suppresses in vitro replication of the R5 strains of HIV-1, which use CCR5 as a coreceptor. Alternative splicing results in multiple transcript variants that encode different isoforms.

Basic Information

Description

Recombinant Human CCL5/RANTES Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Ser24-Ser91) of human CCL5/RANTES (Accession #P13501).

Bio-Activity

Storage

Store the lyophilized protein at -20°C to -80 °C for long term.


After reconstitution, the protein solution is stable at -20 °C for 3 months, at 2-8 °C for up to 1 week.

Avoid repeated freeze/thaw cycles.

Contact

| www.abclonal.com

Validation Data



Recombinant protein Human CCL5/RANTES was determined by SDS-PAGE under reducing conditions with Coomassie Blue, showing a band at 8 kDa.